

# Javani™ Deep Hydrating Gel

## Product Offerings:

- 1 oz. tube (KCJ-501)
- 10-1oz. tubes (KCJ-510)

## Ingredients:

Purified Water, Natural Micronized Collagen, Hyaluronic Acid, Elastin, Sorbitol, Lecithin, Glycerin, Palmitate (Vitamin A), Tocopheryl Acetate (Vitamin E), and Ascorbic Acid (Vitamin C).



## Product Overview:

Javani™ Deep Hydrating Gel is the first hydrating product to deliver micronized, intact collagen that absorbs into the skin. The collagen remains in the skin, drawing in and holding moisture. Javani™ Gel is an excellent choice for deep hydration in a comprehensive skin care program. Javani™ Gel loads skin with moisture, diminishing fine lines and wrinkles and leaving skin with a smooth, supple texture.

- all natural
- alcohol free
- leaves no oily residue
- micronized intact collagen

## Uses:

Javani™ Deep Hydrating Gel can be used where maximum hydration of the skin is necessary. Javani™ Gel can be used in any skin rejuvenation treatment where moisture is lost from the process.

Treatment	Effect
Javani™ Facial	This two-stage facial using Javani™ Gel and Cream makes fine lines disappear, improves texture, and adds a soft, healthy glow.
Laser skin resurfacing	Pre-hydrating the treated areas assists the removal process and prevents excessive moisture loss caused by laser energy.
Microdermabrasion	Hydrating before or after treatment prevents dryness and roughness caused by the removal process.
Chemical peels (includes alpha hydroxy)	Hydrating before treatment prevents moisture loss from deeper layers of the skin. Hydrating after treatment replenishes moisture lost to chemicals used to peel away dead skin cells.

## Application Notes:

Javani™ Deep Hydrating Gel is intended for use by skin care professionals. Javani™ Gel should be applied after cleansing and exfoliation of the area being hydrated. Massage by hand or infuse with an ultrasound or optical machine until completely absorbed. For best results, apply a light coat of Javani™ Gel after it will not absorb any further and cover with a non-absorbing mask. Avoid applying product near the eye and eyelids.

Some larger collagen fibers will not absorb into the skin. Remove these with a damp cloth. Javani™ Gel is not intended as a finishing treatment. It is best used in conjunction with Javani™ Deep Moisturizing Cream.

## Notices:

Contains bovine collagen. Do not use product if you are allergic to bovine materials (e.g. beef products, leather). For external use only. Avoid eye contact. If eye contact occurs, rinse eyes thoroughly with water. Discontinue use if signs of irritation or rash appear. If irritation or rash persists, consult a doctor. Keep out of reach of children.

## Storage:

Store refrigerated between 35-45°F.